

NEW!



E-MAG Active

Mobility with improved safety

ottobock.

Quality for life

E-MAG Active

Mobility with improved safety

The E-MAG Active is used for paralyses of the legs to allow for a physiological gait. The innovative sensor technology, which is also used in smartphones and airplanes, measures the leg's position while walking and controls the orthotic joint accordingly.

New:

Doubly safe: the new E-MAG Active features a PreLock function. This securely locks your knee joint for the stance phase even if you have not fully extended your leg yet. Ask your orthopaedic technician!

Otto Bock HealthCare GmbH · Max-Näder-Straße 15 · 37115 Duderstadt/Germany
T +49 5527 848-1706 · F +49 5527 848-72330 · export@ottobock.de · www.ottobock.com

