


Product Overview

Activation			Stabilisation					
Dorso Carezza Posture	Dorso Direxa Posture		low			medium	high	
50R49	50R59		CI-VII	Lumbo Therna 50R8	Lumbo Carezza 50R40	Lumbo Carezza High 50R42	Lumbo Direxa / Lumbo Direxa Women 50R50 / 50R51	Lumbo Direxa High 50R52
		Thi-XII						
		Lumbar spine L1-V						
+ 2 coil springs	+ 2 plastic rods + 2 metal rods + pad		+ dorsal height: 24 cm + 4 plastic rods	+ dorsal height: 24 cm + 4 plastic rods + optional pad	+ dorsal height: 32 cm + 4 plastic rods + optional pad	+ dorsal height: 24 cm + 4 plastic rods + optional pad + elastic straps + available in waisted version	+ dorsal height: 32 cm + 4 plastic rods + optional pad + elastic straps	+ dorsal height: 28 cm + 4 adaptable metal rods + 2 coil springs + pad + elastic straps



Intelligent, not complicated:

the new Ottobock back orthoses

Ottobock, one of the most innovative healthcare companies in Germany, is now offering a line of back orthoses for all relevant indications. With everything you need for stabilisation and activation in case of back problems.

Enjoy the quick and easy handling of the new back orthoses. See how well your patients feel thanks to the anatomical fit, pleasant wearer comfort and attractive design. We want to make things easier for you and your patients!

Making things easier

► Uniform measurement system

Guideline is the body circumference
2 cm above the iliac crest.

► Uniform sizing system

There are six sizes, each with 10 cm play.
XS is the smallest for a circumference of
70-80 cm, XXL the largest with 120-130 cm.

► Uniform closure system

Practical hand loops for easier application.

Additional advantages

- 1 Textile:** extremely lightweight, easy to clean, machine wash at 40° C
- 2 Design:** modern, sporty
- 3 Pad:** breathable, massaging, stimulates circulation

Quality for life