

ottobock.

Quality for life

5Z15 fibre preform



5Z15 fibre preform

Benefits at a glance

- Quick, clean and easy reinforcement of the socket adapter
- Preform with optimal fibre orientation in the direction of force
- Load-bearing fibre strands stay in place
- Saves time during the reinforcement process
- Very low material waste
- Good draping characteristics
- Simplified quantity control in comparison to rolled goods
- Good absorption of Ottobock lamination resins
- Can be combined with all Ottobock reinforcement materials



The 5Z15 fibre preform reinforcement material was developed especially for the production of fibre-reinforced components using lamination technology. 5Z15 fibre preform is the optimal solution for professionally reinforcing the socket adapter.

Thanks to its specially positioned fibres, the material is designed precisely for the forces at work here. The fibre preform's shape as well as its load-bearing fibre strands increase the efficiency of your workflow.

The lamination process is carried out as usual. Thanks to the fibre preform material, you no longer need to cut reinforcement materials.

Article number	5Z15
Packaging format	2-fibre preform
Order by	Package (2 x fibre preform)
Recommended components*	4R41, 4R42, 4R42=1, 4R43, 4R89, 4R1114, R111=T 4R111=N, 4R116, 4R116=T 4R117, 4R117=T, 4R119=NT, 4R119=T, 4R119, 4R119=N, 4WR95=1, 4WR95=2

*Please note that the 5Z15 fibre preform is incompatible with our socket adapter 5R2, 5R2=C, 4R63, 4R68, 4R100